



# Sterilin Indicator Autoclave Bags

For added assurance Sterilin now offer a high temperature bag with an autoclave indicator. The word 'AUTOCLAVED' appears above 'Hazardous Waste' following a routine autoclave cycle.



# Sterilin Indicator Autoclave Bags



## An advance in laboratory waste disposal

Laboratory waste, which may consist of hazardous chemicals or microorganisms, comes under the Control of Substances Hazardous to Health (COSHH) regulation. This requires that the disposal of such waste is assessed, controlled and suitable precautions put in place to provide adequate protection to the end user.

To support this Sterilin has developed an autoclave bag that has an in-built indicator which clearly tells you that the waste has been autoclaved.

## Features and Benefits

- Unique indicator strip which displays the word 'AUTOCLAVED' following a complete cycle - gives clear indication of successfully autoclaved waste
- Strong disposal bag with blue biohazard printing - ideal for use in all laboratories for hazardous waste disposal
- For use at high temperatures up to 135°C - for the inactivation of particularly resistant waste
- Convenient 'tissue-box' cartons - aids dispensing of the bags
- No need to affix autoclave tape - no training required for new starters in the laboratory.

From a health and safety perspective it is important that all autoclaved waste is clearly identified and segregated; the incorporation of a system that clearly indicates the bag has been autoclaved will fully support biohazard disposal within organisations' waste disposal programs.

## Performance

We can see in the 'before' and 'after' pictures that the integral indicator becomes clearly visible following a routine sterilisation cycle.



## Ordering Information:

Item Code	Description	Max Temp (°C)	Width x Length (mm)	Material	Case Qty
511HTC	Autoclave Bag	135	610 x 810	PP	200



Sterilin is part of Thermo Fisher Scientific

Sterilin Ltd, Parkway, Pen-y-Fan, Newport, NP11 3EF.  
Tel +44 (0)844 844 3737 Fax +44 (0)1495 242 242  
e-mail:info@sterilin.co.uk www.sterilin.co.uk